

产品名称: CDKAL1 Rabbit Polyclonal Antibody
产品货号: APRab08575



产品概述 (Summary)

产品名称 (Production Name)	CDKAL1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDKAL1
别名 (Alternative Names)	CDKAL1; Threonylcarbamoyladenosine tRNA methylthiotransferase; CDK5 regulatory subunit-associated protein 1-like 1; tRNA-t(6)A37 methylthiotransferase
基因 ID (Gene ID)	54901.0
蛋白 ID (SwissProt ID)	Q5VV42.The antiserum was produced against synthesized peptide derived from human CDKAL. AA range:10-59

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	65kDa

产品名称: CDKAL1 Rabbit Polyclonal Antibody
产品货号: APRab08575




研究背景 (Background)

The protein encoded by this gene is a member of the methylthiotransferase family. The function of this gene is not known. Genome-wide association studies have linked single nucleotide polymorphisms in an intron of this gene with susceptibility to type 2 diabetes. [provided by RefSeq, May 2010],cofactor: Binds 2 4Fe-4S clusters. One cluster is coordinated with 3 cysteines and an exchangeable S-adenosyl-L-methionine.,similarity: Belongs to the methylthiotransferase family. CDKAL1 subfamily.,similarity: Contains 1 MTTase N-terminal domain.,similarity: Contains 1 TRAM domain.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of CoLo cells using CDKAL1 Polyclonal Antibody

注意事项 (Note)

For research use only .