

产品名称: CDR1 Rabbit Polyclonal Antibody
产品货号: APRab08584



产品概述 (Summary)

产品名称 (Production Name)	CDR1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDR1
别名 (Alternative Names)	
基因 ID (Gene ID)	1038.0
蛋白 ID (SwissProt ID)	P51861.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	28kDa

研究背景 (Background)

cerebellar degeneration related protein 1(CDR1) Homo sapiens Autoantibodies directed against the protein encoded by this intronless gene have been found in some patients with paraneoplastic cerebellar degeneration. The encoded protein


产品名称: CDR1 Rabbit Polyclonal Antibody
产品货号: APRab08584



contains several hexapeptide repeats. [provided by RefSeq, Jan 2010],disease:CDR1 is one of the target molecules recognized by autoantibodies in patients with paraneoplastic cerebellar degeneration.,sequence caution:Frameshifts in positions 134, 176.,tissue specificity:Brain; predominantly expressed in normal neuroectodermal tissues and in certain malignant tumors.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from KB cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .