

产品名称: CDYL2 Rabbit Polyclonal Antibody
产品货号: APRab08600



产品概述 (Summary)

产品名称 (Production Name)	CDYL2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CDYL2
别名 (Alternative Names)	CDYL2; Chromodomain Y-like protein 2; CDY-like 2
基因 ID (Gene ID)	124359.0
蛋白 ID (SwissProt ID)	Q8N8U2.The antiserum was produced against synthesized peptide derived from human CDYL2. AA range:31-80

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	56kDa

研究背景 (Background)


产品名称: CDYL2 Rabbit Polyclonal Antibody
产品货号: APRab08600




similarity:Contains 1 chromo domain.,tissue specificity:Ubiquitously expressed.,similarity:Contains 1 chromo domain.,tissue specificity:Ubiquitously expressed.,

研究领域 (Research Area)


图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using CDYL2 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using CDYL2 Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .



产品名称: CDYL2 Rabbit Polyclonal Antibody
产品货号: APRab08600



Western Blot analysis of HEPG2 cells using CDYL2 Polyclonal Antibody diluted at 1: 1000 cells nucleus extracted by Minute
TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .