

产品名称: CEACAM21 Rabbit Polyclonal Antibody
产品货号: APRab08616



产品概述 (Summary)

产品名称 (Production Name)	CEACAM21 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEACAM21
别名 (Alternative Names)	CEACAM21; Carcinoembryonic antigen-related cell adhesion molecule 21
基因 ID (Gene ID)	90273.0
蛋白 ID (SwissProt ID)	Q3KPI0.The antiserum was produced against synthesized peptide derived from human CEACAM21. AA range:171-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	32kDa

研究背景 (Background)


产品名称: CEACAM21 Rabbit Polyclonal Antibody
产品货号: APRab08616



similarity:Belongs to the immunoglobulin superfamily. CEA family.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,similarity:Belongs to the immunoglobulin superfamily. CEA family.,similarity:Contains 1 Ig-like C2-type (immunoglobulin-like) domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from HT-29 cells using CEACAM21 antibody.

注意事项 (Note)

For research use only .