

产品名称: CEP152 Rabbit Polyclonal Antibody
产品货号: APRab08659



产品概述 (Summary)

产品名称 (Production Name)	CEP152 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEP152
别名 (Alternative Names)	CEP152; KIAA0912; Centrosomal protein of 152 kDa; Cep152
基因 ID (Gene ID)	22995.0
蛋白 ID (SwissProt ID)	O94986.The antiserum was produced against synthesized peptide derived from human CEP152. AA range:841-890

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	140kDa

研究背景 (Background)


产品名称: CEP152 Rabbit Polyclonal Antibody
产品货号: AP Rab08659




centrosomal protein 152(CEP152) Homo sapiens This gene encodes a protein that is thought to be involved with centrosome function. Mutations in this gene have been associated with primary microcephaly (MCPH4). Alternative splicing results in multiple transcript variants. [provided by RefSeq, Aug 2010],sequence caution:Contaminating sequence. Potential poly-A sequence.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human lymph node tissue, using CEP152 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of HuvEc cells using CEP152 Polyclonal Antibody

注意事项 (Note)

For research use only .