

产品名称: CEP70 Rabbit Polyclonal Antibody
产品货号: APRab08671



产品概述 (Summary)

产品名称 (Production Name)	CEP70 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CEP70
别名 (Alternative Names)	CEP70; BITE; Centrosomal protein of 70 kDa; Cep70; p10-binding protein
基因 ID (Gene ID)	80321.0
蛋白 ID (SwissProt ID)	Q8NHQ1.The antiserum was produced against synthesized peptide derived from human CEP70. AA range:241-290

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	70kDa

研究背景 (Background)


产品名称: CEP70 Rabbit Polyclonal Antibody
产品货号: APRab08671



similarity:Contains 1 TPR repeat.,similarity:Contains 1 TPR repeat.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COS cells, using CEP70 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using CEP70 Polyclonal Antibody diluted at 1: 2000



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtained from antibody

产品名称: CEP70 Rabbit Polyclonal Antibody
产品货号: APRab08671



was pre-absorbed by immunogen peptide.

注意事项 (Note)

For research use only .