

**产品名称: Cerebellin 2 Rabbit Polyclonal Antibody**  
**产品货号: AP Rab08680**



## 产品概述 (Summary)

|                        |                                         |
|------------------------|-----------------------------------------|
| 产品名称 (Production Name) | Cerebellin 2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody              |
| 宿主 (Host)              | Rabbit                                  |
| 应用 (Application)       | WB,ELISA                                |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat,Monkey                  |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                          |
|------------------------|----------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CBLN2                                                                                                    |
| 别名 (Alternative Names) | CBLN2; Cerebellin-2                                                                                      |
| 基因 ID (Gene ID)        | 147381.0                                                                                                 |
| 蛋白 ID (SwissProt ID)   | Q8IUK8.The antiserum was produced against synthesized peptide derived from human CBLN2. AA range:161-210 |

## 产品应用 (Application)

|                          |                                       |
|--------------------------|---------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:20000-1:40000 |
| 蛋白分子量 (Molecular Weight) | 24kDa                                 |

## 研究背景 (Background)


产品名称: Cerebellin 2 Rabbit Polyclonal Antibody  
产品货号: APRab08680




similarity:Contains 1 C1q domain.,subunit:May interact with CBLN1, CBLN3 and CBLN4.,similarity:Contains 1 C1q domain.,subunit:May interact with CBLN1, CBLN3 and CBLN4.,

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of lysates from COS cells, using CBLN2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Cerebellin 2 Polyclonal Antibody

## 注意事项 (Note)

For research use only .