

产品名称: CGI-43 Rabbit Polyclonal Antibody
产品货号: APRab08721



产品概述 (Summary)

产品名称 (Production Name)	CGI-43 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CCZ1/C7orf28B
别名 (Alternative Names)	CCZ1; C7orf28A; CGI-43; Vacuolar fusion protein CCZ1 homolog; CCZ1B; C7orf28B; Vacuolar fusion protein CCZ1 homolog B; Vacuolar fusion protein CCZ1 homolog-like
基因 ID (Gene ID)	221960.0
蛋白 ID (SwissProt ID)	P86791/P86790.The antiserum was produced against synthesized peptide derived from human CG028. AA range:411-460

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	56kDa

产品名称: CGI-43 Rabbit Polyclonal Antibody
产品货号: APRab08721




研究背景 (Background)


similarity:Belongs to the UPF0550 family.,similarity:Belongs to the UPF0550 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat and COS cells, using CG028 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of Jurkat cells using CGI-43 Polyclonal Antibody diluted at 1: 500

注意事项 (Note)

For research use only .