

产品名称: CGR19 Rabbit Polyclonal Antibody
产品货号: APRab08725



产品概述 (Summary)

产品名称 (Production Name)	CGR19 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CGRRF1
别名 (Alternative Names)	CGRRF1; CGR19; RNF197; Cell growth regulator with RING finger domain protein 1; Cell growth regulatory gene 19 protein; RING finger protein 197
基因 ID (Gene ID)	10668.0
蛋白 ID (SwissProt ID)	Q99675.The antiserum was produced against synthesized peptide derived from the C-terminal region of human CGRRF1. AA range:283-332

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	38kDa

产品名称: CGR19 Rabbit Polyclonal Antibody
产品货号: APRab08725




研究背景 (Background)

function:Able to inhibit growth in several cell lines.,similarity:Contains 1 RING-type zinc finger.,function:Able to inhibit growth in several cell lines.,similarity:Contains 1 RING-type zinc finger.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysate from A549 cells, using CGRRF1 Antibody.



Western Blot analysis of A549 cells using CGR19 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .