

产品名称: CGREF1 Rabbit Polyclonal Antibody
产品货号: APRab08727



产品概述 (Summary)

产品名称 (Production Name)	CGREF1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CGREF1
别名 (Alternative Names)	CGREF1; CGR11; Cell growth regulator with EF hand domain protein 1; Cell growth regulatory gene 11 protein; Hydrophobestin
基因 ID (Gene ID)	10669.0
蛋白 ID (SwissProt ID)	Q99674.The antiserum was produced against synthesized peptide derived from the N-terminal region of human CGREF1. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	32kDa

产品名称: CGREF1 Rabbit Polyclonal Antibody
产品货号: APRab08727




研究背景 (Background)


caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Able to inhibit growth in several cell lines.,induction:By p53.,similarity:Contains 2 EF-hand domains.,caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Able to inhibit growth in several cell lines.,induction:By p53.,similarity:Contains 2 EF-hand domains.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysate from 293 cells, using CGREF1 Antibody.




Western Blot analysis of 293 cells using CGREF1 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

产品名称: CGREF1 Rabbit Polyclonal Antibody
产品货号: APRab08727



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100



Immunohistochemical analysis of paraffin-embedded human-brain, antibody was diluted at 1:100

注意事项 (Note)

For research use only .