

产品名称: CHID1 Rabbit Polyclonal Antibody
产品货号: APRab08759



产品概述 (Summary)

产品名称 (Production Name)	CHID1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHID1 GL008 PSEC0104 SB139
别名 (Alternative Names)	CHID1 GL008 PSEC0104 SB139
基因 ID (Gene ID)	66005.0
蛋白 ID (SwissProt ID)	Q9BWS9.The antiserum was produced against synthesized peptide derived from the Internal region of human CHID1. AA range:81-130

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	45kDa

研究背景 (Background)


产品名称: CHID1 Rabbit Polyclonal Antibody
产品货号: APRab08759



induction:Up-regulated by interleukin-4 and dexamethasone in the macrophages.,similarity:Belongs to the glycosyl hydrolase 18 family.,subunit:Interacts with STAB1.,tissue specificity:Expressed in cells of monocytic, T, B and epithelial origin.,induction:Up-regulated by interleukin-4 and dexamethasone in the macrophages.,similarity:Belongs to the glycosyl hydrolase 18 family.,subunit:Interacts with STAB1.,tissue specificity:Expressed in cells of monocytic, T, B and epithelial origin.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of 3T3 cells using Antibody diluted at 500. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .