

产品名称: CHSTE Rabbit Polyclonal Antibody
产品货号: APRab08795



产品概述 (Summary)

产品名称 (Production Name)	CHSTE Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CHST14 D4ST1 UNQ1925/PRO4400
别名 (Alternative Names)	
基因 ID (Gene ID)	113189.0
蛋白 ID (SwissProt ID)	Q8NCH0.Synthesized peptide derived from part region of human protein

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	41kDa

研究背景 (Background)

This gene encodes a member of the HNK-1 family of sulfotransferases. The encoded protein transfers sulfate to the C-4 hydroxyl of N-acetylgalactosamine residues in dermatan sulfate. Mutations in this gene have been associated with

产品名称: CHSTE Rabbit Polyclonal Antibody
产品货号: APRab08795




adducted thumb-clubfoot syndrome.[provided by RefSeq, Mar 2010],function:Catalyzes the transfer of sulfate to position 4 of the N-acetylgalactosamine (GalNAc) residue of dermatan sulfate. Transfers sulfate to the C-4 hydroxyl of beta1,4-linked GalNAc that is substituted with an alpha-linked iduronic acid (IdoUA) at the C-3 hydroxyl. Transfers sulfate more efficiently to GalNAc residues in -IdoUA-GalNAc-IdoUA- than in -GlcUA-GalNAc-GlcUA-sequences. Has preference for partially desulfated dermatan sulfate. Addition of sulfate to GalNAc may occur immediately after epimerization of GlcUA to IdoUA.,similarity:Belongs to the sulfotransferase 2 family.,tissue specificity:Widely expressed. Expressed at high level in pituitary gland, placenta, uterus and thyroid.,

研究领域 (Research Area)

Chondroitin sulfate biosynthesis;

图片 (Image Data)



Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .