

产品名称: CIB3 Rabbit Polyclonal Antibody
产品货号: APRab08809



产品概述 (Summary)

产品名称 (Production Name)	CIB3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CIB3
别名 (Alternative Names)	CIB3; KIP3; Calcium and integrin-binding family member 3; Kinase-interacting protein 3; KIP 3
基因 ID (Gene ID)	117286.0
蛋白 ID (SwissProt ID)	Q96Q77.The antiserum was produced against synthesized peptide derived from human CIB3. AA range:138-187

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	32kDa

产品名称: CIB3 Rabbit Polyclonal Antibody
产品货号: APRab08809




研究背景 (Background)


This gene product shares a high degree of sequence similarity with DNA-dependent protein kinase catalytic subunit-interacting protein 2 in human and mouse, and like them may bind the catalytic subunit of DNA-dependent protein kinases. The exact function of this gene is not known. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2014], similarity: Contains 3 EF-hand domains.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from K562 cells, using CIB3 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of K562 cells using CIB3 Polyclonal Antibody

注意事项 (Note)

For research use only .