

产品名称: CIDE-B Rabbit Polyclonal Antibody
产品货号: APRab08814



产品概述 (Summary)

产品名称 (Production Name)	CIDE-B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CIDEB
别名 (Alternative Names)	CIDEB; Cell death activator CIDE-B; Cell death-inducing DFFA-like effector B
基因 ID (Gene ID)	27141.0
蛋白 ID (SwissProt ID)	Q9UHD4.The antiserum was produced against synthesized peptide derived from human CIDEB. AA range:91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)


产品名称: CIDE-B Rabbit Polyclonal Antibody
产品货号: APRab08814




function:Activates apoptosis.,similarity:Contains 1 CIDE-N domain.,subunit:Inhibited by DFFB. Interacts with DFFA and DFFB.,tissue specificity:Highly expressed in liver and small intestine and, at lower levels, in colon, kidney and spleen.,function:Activates apoptosis.,similarity:Contains 1 CIDE-N domain.,subunit:Inhibited by DFFB. Interacts with DFFA and DFFB.,tissue specificity:Highly expressed in liver and small intestine and, at lower levels, in colon, kidney and spleen.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of COS7 cells, using CIDEB Antibody. The picture on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human cervix caecinoma tissue, using CIDEB Antibody. The picture on the right is blocked with the synthesized peptide.

产品名称: CIDE-B Rabbit Polyclonal Antibody
产品货号: APRab08814



Western blot analysis of lysates from HT-29 and RAW264.7 cells, using CIDE-B Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using CIDE-B Polyclonal Antibody

注意事项 (Note)

For research use only .