

产品名称: CLCC1 Rabbit Polyclonal Antibody
产品货号: APRab08930



产品概述 (Summary)

产品名称 (Production Name)	CLCC1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLCC1
别名 (Alternative Names)	CLCC1; KIAA0761; MCLC; Chloride channel CLIC-like protein 1; Mid-1-related chloride channel protein 1
基因 ID (Gene ID)	23155.0
蛋白 ID (SwissProt ID)	Q96S66. The antiserum was produced against synthesized peptide derived from human CLCC1. AA range: 391-440

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:200-1:1000, ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	62kDa

产品名称: CLCC1 Rabbit Polyclonal Antibody
产品货号: APRab08930




研究背景 (Background)


function: Seems to act as a chloride ion channel., similarity: Belongs to the chloride channel MCLC family., function: Seems to act as a chloride ion channel., similarity: Belongs to the chloride channel MCLC family.,

研究领域 (Research Area)

图片 (Image Data)




Immunofluorescence analysis of A549 cells, using CLCC1 Antibody. The picture on the right is blocked with the synthesized




Western blot analysis of lysates from COLO and HepG2 cells, using CLCC1 Antibody. The lane on the right is blocked with the synthesized peptide.


产品名称: CLCC1 Rabbit Polyclonal Antibody
产品货号: APRab08930



Western blot analysis of the lysates from HT-29 cells using CLCC1 antibody.



Western Blot analysis of various cells using CLCC1 Polyclonal Antibody diluted at 1: 1000



Western Blot analysis of HepG2 cells using CLCC1 Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .