

**产品名称: Cleaved-CD97 $\alpha$  (L530) Rabbit Polyclonal Antibody**  
**产品货号: APRab08980**



## 产品概述 (Summary)

|                        |                                                         |
|------------------------|---------------------------------------------------------|
| 产品名称 (Production Name) | Cleaved-CD97 $\alpha$ (L530) Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody                              |
| 宿主 (Host)              | Rabbit                                                  |
| 应用 (Application)       | WB,ELISA                                                |
| 种属反应性 (Reactivity)     | Human,Rat,Mouse                                         |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                              |
|------------------------|--------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CD97                                                                                                         |
| 别名 (Alternative Names) | CD97; CD97 antigen; Leukocyte antigen CD97; CD antigen CD97                                                  |
| 基因 ID (Gene ID)        | 976.0                                                                                                        |
| 蛋白 ID (SwissProt ID)   | P48960.The antiserum was produced against synthesized peptide derived from human CD97alpha. AA range:481-530 |

## 产品应用 (Application)

|                          |                                      |
|--------------------------|--------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,ELISA 1:5000-1:10000 |
| 蛋白分子量 (Molecular Weight) | 55kDa                                |


## 研究背景 (Background)

**产品名称: Cleaved-CD97 $\alpha$  (L530) Rabbit Polyclonal Antibody**  
**产品货号: APRab08980**

This gene encodes a member of the EGF-TM7 subfamily of adhesion G protein-coupled receptors, which mediate cell-cell interactions. These proteins are cleaved by self-catalytic proteolysis into a large extracellular subunit and seven-span transmembrane subunit, which associate at the cell surface as a receptor complex. The encoded protein may play a role in cell adhesion as well as leukocyte recruitment, activation and migration, and contains multiple extracellular EGF-like repeats which mediate binding to chondroitin sulfate and the cell surface complement regulatory protein CD55. Expression of this gene may play a role in the progression of several types of cancer. Alternatively spliced transcript variants encoding multiple isoforms with 3 to 5 EGF-like repeats have been observed for this gene. This gene is found in a cluster with other EGF-TM7 genes on the short arm of chromosome 1. domain: Binding to chondroitin sulfate is mediated by the fourth EGF domain. domain: The first two EGF domains mediate the interaction with DAF. A third tandemly arranged EGF domain is necessary for the structural integrity of the binding region. function: Receptor potentially involved in both adhesion and signaling processes early after leukocyte activation. Plays an essential role in leukocyte migration. induction: Rapid up-regulation during lymphocyte activation. PTM: Proteolytically cleaved into 2 subunits, an extracellular alpha subunit and a seven-transmembrane subunit. similarity: Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily. similarity: Contains 1 GPS domain. similarity: Contains 5 EGF-like domains. subunit: Forms a heterodimer, consisting of a large extracellular region (alpha subunit) non-covalently linked to a seven-transmembrane moiety (beta subunit). Interacts with complement decay-accelerating factor (DAF). The largest isoform (isoform 1) interacts with chondroitin sulfate. tissue specificity: Broadly expressed, found on most hematopoietic cells, including activated lymphocytes, monocytes, macrophages, dendritic cells, and granulocytes. Expressed also abundantly by smooth muscle cells. Expressed in thyroid, colorectal, gastric, esophageal and pancreatic carcinomas too. Expression are increased under inflammatory conditions in the CNS of multiple sclerosis and in synovial tissue of patients with rheumatoid arthritis. Increased expression of CD97 in the synovium is accompanied by detectable levels of soluble CD97 in the synovial fluid.

## 研究领域 (Research Area)

## 图片 (Image Data)



Western blot analysis of CD97 $\alpha$  (Cleaved-Leu530) Antibody. The lane on the right is blocked with the CD97 $\alpha$  antibody.

**产品名称: Cleaved-CD97 $\alpha$  (L530) Rabbit Polyclonal Antibody**  
**产品货号: APRab08980**



---

(Cleaved-Leu530) peptide.

**注意事项 (Note)**

For research use only .