

产品名称: Cleaved-Ephrin-A2 (N188) Rabbit Polyclonal Antibody
产品货号: APRab08985



产品概述 (Summary)

产品名称 (Production Name)	Cleaved-Ephrin-A2 (N188) Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EFNA2
别名 (Alternative Names)	EFNA2; EPLG6; LERK6; Ephrin-A2; EPH-related receptor tyrosine kinase ligand 6; LERK-6; HEK7 ligand; HEK7-L
基因 ID (Gene ID)	1943.0
蛋白 ID (SwissProt ID)	O43921.The antiserum was produced against synthesized peptide derived from human EFNA2. AA range:139-188

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	21kDa

产品名称: Cleaved-Ephrin-A2 (N188) Rabbit Polyclonal Antibody
产品货号: AP Rab08985




研究背景 (Background)

This gene encodes a member of the ephrin family. The protein is composed of a signal sequence, a receptor-binding region, a spacer region, and a hydrophobic region. The EPH and EPH-related receptors comprise the largest subfamily of receptor protein-tyrosine kinases and have been implicated in mediating developmental events, particularly in the nervous system. Based on their structures and sequence relationships, ephrins are divided into the ephrin-A (EFNA) class, which are anchored to the membrane by a glycosylphosphatidylinositol linkage, and the ephrin-B (EFNB) class, which are transmembrane proteins. Posttranslational modifications determine whether this protein localizes to the nucleus or the cytoplasm. [provided by RefSeq, Jul 2008], similarity: Belongs to the ephrin family., subunit: Binds to the receptor tyrosine kinases EPHA3, EPHA4 and EPHA5.,

研究领域 (Research Area)


Axon guidance;

图片 (Image Data)



Western blot analysis of lysates from 293 cells, treated with etoposide 25uM 1h, using EFNA2 (Cleaved-Asn188) Antibody.

The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using Cleaved-Ephrin-A2 (N188) Polyclonal Antibody

产品名称: Cleaved-Ephrin-A2 (N188) Rabbit Polyclonal
Antibody
产品货号: APRab08985



注意事项 (Note)

For research use only .