

产品名称: CLIP-170 Rabbit Polyclonal Antibody
产品货号: APRab09045



产品概述 (Summary)

产品名称 (Production Name)	CLIP-170 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLIP1 CLIP1; CYLN1; RSN; CAP-Gly domain-containing linker protein 1; Cytoplasmic linker protein 1; Cytoplasmic linker protein 170 alpha-2; CLIP-170; Reed-Sternberg intermediate filament-associated protein; Restin
别名 (Alternative Names)	
基因 ID (Gene ID)	6249.0
蛋白 ID (SwissProt ID)	P30622.The antiserum was produced against synthesized peptide derived from human CLIP1. AA range:1291-1340

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	161kDa

产品名称: CLIP-170 Rabbit Polyclonal Antibody
产品货号: APRab09045




研究背景 (Background)

The protein encoded by this gene links endocytic vesicles to microtubules. This gene is highly expressed in Reed-Sternberg cells of Hodgkin disease. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],function:Seems to be a intermediate filament associated protein that links endocytic vesicles to microtubules.,similarity:Contains 2 CAP-Gly domains.,subcellular location:Associated with the cytoskeleton.,tissue specificity:Highly expressed in the Reed-Sternberg cells of Hodgkin's disease.,


研究领域 (Research Area)

Regulation of Microtubule Dynamics

图片 (Image Data)




Western blot analysis of lysates from HUVEC, COLO, MCF-7, HepG2, and RAW264.7 cells, using CLIP1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using CLIP-170 Polyclonal Antibody

产品名称: CLIP-170 Rabbit Polyclonal Antibody
产品货号: APRab09045



Western Blot analysis of RAW264.7 cells using CLIP-170 Polyclonal Antibody

注意事项 (Note)

For research use only .