

产品名称: CLIP4 Rabbit Polyclonal Antibody
产品货号: APRab09047



产品概述 (Summary)

产品名称 (Production Name)	CLIP4 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CLIP4
别名 (Alternative Names)	CLIP4; RSNL2; CAP-Gly domain-containing linker protein 4; Restin-like protein 2
基因 ID (Gene ID)	79745.0
蛋白 ID (SwissProt ID)	Q8N3C7.The antiserum was produced against synthesized peptide derived from human CLIP4. AA range:321-370

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	76kDa

产品名称: CLIP4 Rabbit Polyclonal Antibody
产品货号: APRab09047




研究背景 (Background)

similarity:Contains 3 ANK repeats.,similarity:Contains 3 CAP-Gly domains.,similarity:Contains 3 ANK repeats.,similarity:Contains 3 CAP-Gly domains.,


研究领域 (Research Area)

Regulation of Microtubule Dynamics


图片 (Image Data)



Western blot analysis of CLIP4 Antibody. The lane on the right is blocked with the CLIP4 peptide.



Western blot analysis of the lysates from HUVEC cells using CLIP4 antibody.



Western Blot analysis of various cells using CLIP4 Polyclonal Antibody

产品名称: CLIP4 Rabbit Polyclonal Antibody
产品货号: APRab09047



注意事项 (Note)

For research use only .