

产品名称: CMTM2 Rabbit Polyclonal Antibody
产品货号: APRab09085



产品概述 (Summary)

产品名称 (Production Name)	CMTM2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CMTM2
别名 (Alternative Names)	CMTM2; CKLFSF2; CKLF-like MARVEL transmembrane domain-containing protein 2; Chemokine-like factor superfamily member 2
基因 ID (Gene ID)	146225.0
蛋白 ID (SwissProt ID)	Q8TAZ6.The antiserum was produced against synthesized peptide derived from human CKLF2. AA range:171-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	27kDa

产品名称: CMTM2 Rabbit Polyclonal Antibody
产品货号: APRab09085




研究背景 (Background)


This gene belongs to the chemokine-like factor gene superfamily, a novel family that links the chemokine and the transmembrane 4 superfamilies of signaling molecules. The protein encoded by this gene may play an important role in testicular development. [provided by RefSeq, Jul 2008],similarity:Belongs to the chemokine-like factor family.,similarity:Contains 1 MARVEL domain.,tissue specificity:Highly expressed in testis.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from K562 cells, using CKLF2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using CKLF2 antibody.

产品名称: CMTM2 Rabbit Polyclonal Antibody
产品货号: APRab09085



Western Blot analysis of K562 cells using CMTM2 Polyclonal Antibody

注意事项 (Note)

For research use only .