

产品名称: CPXM2 Rabbit Polyclonal Antibody  
产品货号: APRab09338



## 产品概述 (Summary)

|                        |                                  |
|------------------------|----------------------------------|
| 产品名称 (Production Name) | CPXM2 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody       |
| 宿主 (Host)              | Rabbit                           |
| 应用 (Application)       | IHC, ICC, IF, ELISA              |
| 种属反应性 (Reactivity)     | Human, Rat, Mouse                |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                            |
|------------------------|------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CPXM2                                                                                                      |
| 别名 (Alternative Names) | CPXM2; Inactive carboxypeptidase-like protein X2                                                           |
| 基因 ID (Gene ID)        | 119587.0                                                                                                   |
| 蛋白 ID (SwissProt ID)   | Q8N436. The antiserum was produced against synthesized peptide derived from human CPXM2. AA range: 450-500 |

## 产品应用 (Application)

|                          |                                                           |
|--------------------------|-----------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) |                                                           |

## 研究背景 (Background)


**产品名称: CPXM2 Rabbit Polyclonal Antibody**  
**产品货号: APRab09338**



caution:As it has lost active sites residues and zinc-binding sites it is unlikely to be catalytically active.,function:May be involved in cell-cell interactions.,similarity:Belongs to the peptidase M14 family.,similarity:Contains 1 F5/8 type C domain.,caution:As it has lost active sites residues and zinc-binding sites it is unlikely to be catalytically active.,function:May be involved in cell-cell interactions.,similarity:Belongs to the peptidase M14 family.,similarity:Contains 1 F5/8 type C domain.,

## **研究领域 (Research Area)**

## **图片 (Image Data)**



Immunohistochemical analysis of paraffin-embedded human tonsil. 1, Antibody was diluted at 1:200 (4° overnight) . 2, Tris-EDTA,pH9.0 was used for antigen retrieval. 3,Secondary antibody was diluted at 1:200 (room temperature, 30min) .

## **注意事项 (Note)**

For research use only .