

产品名称: CREM Rabbit Polyclonal Antibody  
产品货号: APRab09388



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | CREM Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA             |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                 |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                        |
|------------------------|--------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | CREM                                                                                                   |
| 别名 (Alternative Names) | CREM; cAMP-responsive element modulator; Inducible cAMP early repressor; ICER                          |
| 基因 ID (Gene ID)        | 1390.0                                                                                                 |
| 蛋白 ID (SwissProt ID)   | Q03060.The antiserum was produced against synthesized peptide derived from human CREM. AA range:81-130 |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 34kDa                                                                   |

产品名称: CREM Rabbit Polyclonal Antibody  
产品货号: APRab09388




## 研究背景 (Background)

cAMP responsive element modulator(CREM) Homo sapiens This gene encodes a bZIP transcription factor that binds to the cAMP responsive element found in many viral and cellular promoters. It is an important component of cAMP-mediated signal transduction during the spermatogenic cycle, as well as other complex processes. Alternative promoter and translation initiation site usage allows this gene to exert spatial and temporal specificity to cAMP responsiveness. Multiple alternatively spliced transcript variants encoding several different isoforms have been found for this gene, with some of them functioning as activators and some as repressors of transcription. [provided by RefSeq, Jul 2008],caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,function:Transcriptional regulator that binds the cAMP response element (CRE), a sequence present in many viral and cellular promoters. Isoforms are either transcriptional activators or repressors. Plays a role in spermatogenesis and is involved in spermatid maturation.,PTM:Stimulated by phosphorylation.,sequence caution:Translated as Tyr.,similarity:Belongs to the bZIP family.,similarity:Contains 1 bZIP domain.,similarity:Contains 1 KID (kinase-inducible) domain.,subunit:Binds DNA as a dimer. Interacts with FHL5.,tissue specificity:Expressed in testes (round spermatids) (at protein level). Isoform 14 is the major activator form in testes.,


## 研究领域 (Research Area)

## 图片 (Image Data)




Western blot analysis of lysates from RAW264.7 cells, using CREM Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: CREM Rabbit Polyclonal Antibody  
产品货号: APRab09388



Immunohistochemical analysis of paraffin-embedded Human brain. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negative contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.



Immunohistochemical analysis of paraffin-embedded Human breast cancer. Antibody was diluted at 1:100 (4°,overnight) . High-pressure and temperature Tris-EDTA,pH8.0 was used for antigen retrieval. Negetive contrl (right) obtaned from antibody was pre-absorbed by immunogen peptide.

### 注意事项 (Note)

For research use only .