

产品名称: CSRNP3 Rabbit Polyclonal Antibody
产品货号: APRab09468



产品概述 (Summary)

产品名称 (Production Name)	CSRNP3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CSRNP3
别名 (Alternative Names)	CSRNP3; FAM130A2; TAIP2; Cysteine/serine-rich nuclear protein 3; CSRNP-3; Protein FAM130A2; TGF-beta-induced apoptosis protein 2; TAIP-2
基因 ID (Gene ID)	80034.0
蛋白 ID (SwissProt ID)	Q8WYN3.The antiserum was produced against synthesized peptide derived from human TAIP-2. AA range:171-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	64kDa

产品名称: CSRNP3 Rabbit Polyclonal Antibody
产品货号: APRab09468




研究背景 (Background)

function: Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis., sequence caution: Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon., similarity: Belongs to the AXUD1 family., function: Binds to the consensus sequence 5'-AGAGTG-3' and has transcriptional activator activity. Plays a role in apoptosis., sequence caution: Unusual initiator. The initiator methionine is coded by a non-canonical CTG leucine codon., similarity: Belongs to the AXUD1 family.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of the lysates from HepG2 cells using TAIP-2 antibody.

注意事项 (Note)

For research use only .