

产品名称: CstF-64T Rabbit Polyclonal Antibody
产品货号: APRab09475



产品概述 (Summary)

产品名称 (Production Name)	CstF-64T Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CSTF2T
别名 (Alternative Names)	CSTF2T; KIAA0689; Cleavage stimulation factor subunit 2 tau variant; CF-1 64 kDa subunit tau variant; Cleavage stimulation factor 64 kDa subunit tau variant; CSTF 64 kDa subunit tau variant; TauCstF-64
基因 ID (Gene ID)	23283.0
蛋白 ID (SwissProt ID)	Q9H0L4.The antiserum was produced against synthesized peptide derived from human CSTF2T. AA range:91-140

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	64kDa

产品名称: CstF-64T Rabbit Polyclonal Antibody
产品货号: APRab09475




研究背景 (Background)


function: May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 1 RRM (RNA recognition motif) domain., function: May play a significant role in AAUAAA-independent mRNA polyadenylation in germ cells. Directly involved in the binding to pre-mRNAs., PTM: Phosphorylated upon DNA damage, probably by ATM or ATR., similarity: Contains 1 RRM (RNA recognition motif) domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from Jurkat and MCF-7 cells, using CSTF2T Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from Jurkat cells using CSTF2T antibody.

产品名称: CstF-64T Rabbit Polyclonal Antibody
产品货号: APRab09475



Western Blot analysis of Jurkat cells using CstF-64T Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventbiotech, MN, USA) .

注意事项 (Note)

For research use only .