

产品名称: CYB561D2 Rabbit Polyclonal Antibody
产品货号: APRab09571



产品概述 (Summary)

产品名称 (Production Name)	CYB561D2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYB561D2
别名 (Alternative Names)	CYB561D2; 101F6; LUCA12.2; Cytochrome b561 domain-containing protein 2; Putative tumor suppressor protein 101F6
基因 ID (Gene ID)	11068.0
蛋白 ID (SwissProt ID)	O14569.The antiserum was produced against synthesized peptide derived from human C56D2. AA range:55-104

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

产品名称: CYB561D2 Rabbit Polyclonal Antibody
产品货号: APRab09571




研究背景 (Background)


cofactor: Binds 2 heme groups non-covalently., similarity: Contains 1 cytochrome b561 domain., cofactor: Binds 2 heme groups non-covalently., similarity: Contains 1 cytochrome b561 domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using C56D2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of Jurkat cells using CYB561D2 Polyclonal Antibody

注意事项 (Note)

For research use only .