

产品名称: CYP2A6V2 Rabbit Polyclonal Antibody
产品货号: APRab09648



产品概述 (Summary)

产品名称 (Production Name)	CYP2A6V2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYP2A6V2
别名 (Alternative Names)	
基因 ID (Gene ID)	1548.0
蛋白 ID (SwissProt ID)	Q13120.Synthesized peptide derived from the Internal region of human CYP2A6V2.

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	56kDa

研究背景 (Background)

产品名称: CYP2A6V2 Rabbit Polyclonal Antibody
产品货号: APRab09648




This gene, CYP2A6, encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. This protein localizes to the endoplasmic reticulum and its expression is induced by phenobarbital. The enzyme is known to hydroxylate coumarin, and also metabolizes nicotine, aflatoxin B1, nitrosamines, and some pharmaceuticals. Individuals with certain allelic variants are said to have a poor metabolizer phenotype, meaning they do not efficiently metabolize coumarin or nicotine. This gene is part of a large cluster of cytochrome P450 genes from the CYP2A, CYP2B and CYP2F subfamilies on chromosome 19q. The gene was formerly referred to as CYP2A3; however, it has been renamed CYP2A6. oxidation reduction,

研究领域 (Research Area)

Caffeine metabolism;Retinol metabolism;Drug metabolism;Drug metabolism;

图片 (Image Data)



Western Blot analysis of L929 cells using CYP2A6V2 Polyclonal Antibody diluted at 1: 2000

注意事项 (Note)

For research use only .