

产品名称: CYTL1 Rabbit Polyclonal Antibody
产品货号: APRab09704



产品概述 (Summary)

产品名称 (Production Name)	CYTL1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	CYTL1
别名 (Alternative Names)	CYTL1; C4orf4; Cytokine-like protein 1; Protein C17
基因 ID (Gene ID)	54360.0
蛋白 ID (SwissProt ID)	Q9NRR1.The antiserum was produced against synthesized peptide derived from human CYTL1. AA range:61-110

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	16kDa

研究背景 (Background)


产品名称: CYTL1 Rabbit Polyclonal Antibody
产品货号: AP Rab09704




C17 is a cytokine-like protein specifically expressed in bone marrow and cord blood mononuclear cells that bear the CD34 (MIM 142230) surface marker (Liu et al., 2000 [PubMed 10857752]).[supplied by OMIM, Mar 2008],tissue specificity:Specifically expressed in CD34+ hematopoietic cells.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human brain tissue, using CYTL1 Antibody. The picture on the right is blocked with the synthesized peptide.



Western blot analysis of lysates from LOVO cells, using CYTL1 Antibody. The lane on the right is blocked with the synthesized peptide.

产品名称: CYTL1 Rabbit Polyclonal Antibody
产品货号: APRab09704



Western blot analysis of the lysates from HT-29 cells using CYTL1 antibody.

注意事项 (Note)

For research use only .