

产品名称: D54 Rabbit Polyclonal Antibody
产品货号: APRab09771



产品概述 (Summary)

产品名称 (Production Name)	D54 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	TPD52L2
别名 (Alternative Names)	TPD52L2; Tumor protein D54; hD54; Tumor protein D52-like 2
基因 ID (Gene ID)	7165.0
蛋白 ID (SwissProt ID)	O43399.The antiserum was produced against synthesized peptide derived from human TPD54. AA range:141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	22kDa

研究背景 (Background)


产品名称: D54 Rabbit Polyclonal Antibody
产品货号: APRab09771




This gene encodes a member of the tumor protein D52-like family. These proteins are characterized by an N-terminal coiled-coil motif that is used to form homo- and heteromeric complexes with other tumor protein D52-like proteins. Expression of this gene may be a marker for breast cancer and acute lymphoblastic leukemia. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene, and a pseudogene of this gene is located on the long arm of chromosome 12. [provided by RefSeq, Aug 2011],alternative products:Additional isoforms seem to exist,similarity:Belongs to the TPD52 family.,subunit:Forms homodimer or heterodimer with other members of the family (By similarity). Interacts with MAL2.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using TPD54 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using D54 Polyclonal Antibody

注意事项 (Note)

For research use only .