

产品名称: DCAMKL3 Rabbit Polyclonal Antibody
产品货号: AP Rab09830



产品概述 (Summary)

产品名称 (Production Name)	DCAMKL3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DCLK3 DCLK3; DCAMKL3; DCDC3C; KIAA1765; Serine/threonine-protein kinase
别名 (Alternative Names)	DCLK3; Doublecortin domain-containing protein 3C; Doublecortin-like and CAM kinase-like 3; Doublecortin-like kinase 3
基因 ID (Gene ID)	85443.0
蛋白 ID (SwissProt ID)	Q9C098.The antiserum was produced against synthesized peptide derived from human DCLK3. AA range:291-340

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	73kDa

产品名称: DCAMKL3 Rabbit Polyclonal Antibody
产品货号: APRab09830




研究背景 (Background)


catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. CaMK subfamily.,similarity:Contains 1 protein kinase domain.,catalytic activity:ATP + a protein = ADP + a phosphoprotein.,similarity:Belongs to the protein kinase superfamily.,similarity:Belongs to the protein kinase superfamily. CAMK Ser/Thr protein kinase family. CaMK subfamily.,similarity:Contains 1 protein kinase domain.,

研究领域 (Research Area)

图片 (Image Data)




Immunohistochemistry analysis of paraffin-embedded human brain, using DCLK3 Antibody. The picture on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using DCAMKL3 Polyclonal Antibody

产品名称: DCAMKL3 Rabbit Polyclonal Antibody
产品货号: APRab09830



Western Blot analysis of Jurkat cells using DCAMKL3 Polyclonal Antibody

注意事项 (Note)

For research use only .