

产品名称: DDT Rabbit Polyclonal Antibody
产品货号: APRab09870



产品概述 (Summary)

产品名称 (Production Name)	DDT Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse, Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DDT
别名 (Alternative Names)	D-dopachrome decarboxylase (EC 4.1.1.84; D-dopachrome tautomerase; Phenylpyruvate tautomerase II)
基因 ID (Gene ID)	1652.0
蛋白 ID (SwissProt ID)	P30046.Synthetic peptide from human protein at AA range: 50-90

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: DDT Rabbit Polyclonal Antibody
产品货号: APRab09870




D-dopachrome tautomerase converts D-dopachrome into 5,6-dihydroxyindole. The DDT gene is related to the migration inhibitory factor (MIF) in terms of sequence, enzyme activity, and gene structure. DDT and MIF are closely linked on chromosome 22. [provided by RefSeq, Jul 2008], catalytic activity: D-dopachrome = 5,6-dihydroxyindole + CO(2), function: Tautomerization of D-dopachrome with decarboxylation to give 5,6-dihydroxyindole (DHI), similarity: Belongs to the MIF family, subunit: Homotrimer,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-placenta, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-placenta, antibody was diluted at 1:200

注意事项 (Note)

For research use only .