

产品名称: Defensin α 3 Rabbit Polyclonal Antibody
产品货号: APRab09903



产品概述 (Summary)

产品名称 (Production Name)	Defensin α 3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DEFA3
别名 (Alternative Names)	DEFA3; DEF3; Neutrophil defensin 3; Defensin; α 3; HNP-3; HP-3; HP3
基因 ID (Gene ID)	1668.0
蛋白 ID (SwissProt ID)	P59666. The antiserum was produced against synthesized peptide derived from human Defensin α 3. AA range: 3-52

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: Defensin α 3 Rabbit Polyclonal Antibody
产品货号: APRab09903



defensin alpha 3(DEF3) Homo sapiens Defensins are a family of antimicrobial and cytotoxic peptides thought to be involved in host defense. They are abundant in the granules of neutrophils and also found in the epithelia of mucosal surfaces such as those of the intestine, respiratory tract, urinary tract, and vagina. Members of the defensin family are highly similar in protein sequence and distinguished by a conserved cysteine motif. The protein encoded by this gene, defensin, alpha 3, is found in the microbicidal granules of neutrophils and likely plays a role in phagocyte-mediated host defense. Several alpha defensin genes are clustered on chromosome 8. This gene differs from defensin, alpha 1 by only one amino acid. This gene and the gene encoding defensin, alpha 1 are both subject to copy number variation. [provided by RefSeq, Oct 2014],function:Defensin 2 and defensin 3 have antibiotic, fungicide and antiviral activities. Has antimicrobial activity against Gram-negative and Gram-positive bacteria. Defensins are thought to kill microbes by permeabilizing their plasma membrane.,similarity:Belongs to the alpha-defensin family.,subunit:Dimer.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of Defensin α 3 antibody in paraffin-embedded human skeletal muscle tissue.

注意事项 (Note)

For research use only .