

产品名称: DGAT2L6 Rabbit Polyclonal Antibody
产品货号: APRab09940



产品概述 (Summary)

产品名称 (Production Name)	DGAT2L6 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DGAT2L6
别名 (Alternative Names)	DGAT2L6; DC3; Diacylglycerol O-acyltransferase 2-like protein 6; Diacylglycerol O-acyltransferase candidate 3; hDC3
基因 ID (Gene ID)	347516.0
蛋白 ID (SwissProt ID)	Q6ZPD8.The antiserum was produced against synthesized peptide derived from human DGAT2L6. AA range:271-320

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	36kDa

产品名称: DGAT2L6 Rabbit Polyclonal Antibody
产品货号: APRab09940




研究背景 (Background)

This gene is a member of the diacylglycerol acyltransferase 2 family. The encoded protein is a putative acyltransferase and is most likely involved in the synthesis of di- or triacylglycerol, however its substrate specificity is currently unknown. [provided by RefSeq, Feb 2015],function:Probable acyltransferase uses fatty acyl-CoA as substrate (By similarity). Has no wax synthase activity to produce wax esters.,similarity:Belongs to the diacylglycerol acyltransferase family.,tissue specificity:Expressed in all tissues tested except pancreas.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from Jurkat cells, using DGAT2L6 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .