

产品名称: DGK- β Rabbit Polyclonal Antibody
产品货号: APRab09947



产品概述 (Summary)

产品名称 (Production Name)	DGK- β Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DGKB
别名 (Alternative Names)	DGKB; DAGK2; KIAA0718; Diacylglycerol kinase beta; DAG kinase beta; 90 kDa diacylglycerol kinase; Diglyceride kinase beta; DGK-beta
基因 ID (Gene ID)	1607.0
蛋白 ID (SwissProt ID)	Q9Y6T7.The antiserum was produced against synthesized peptide derived from human DGKB. AA range:657-706

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	90kDa

产品名称: DGK- β Rabbit Polyclonal Antibody
产品货号: APRab09947




研究背景 (Background)

Diacylglycerol kinases (DGKs) are regulators of the intracellular concentration of the second messenger diacylglycerol (DAG) and thus play a key role in cellular processes. Nine mammalian isotypes have been identified, which are encoded by separate genes. Mammalian DGK isozymes contain a conserved catalytic (kinase) domain and a cysteine-rich domain (CRD). The protein encoded by this gene is a diacylglycerol kinase, beta isotype. Two alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008],catalytic activity:ATP + 1,2-diacylglycerol = ADP + 1,2-diacyl-sn-glycerol 3-phosphate.,enzyme regulation:Stimulated by phosphatidylserine.,function:Exhibits high phosphorylation activity for long-chain diacylglycerols.,similarity:Belongs to the eukaryotic diacylglycerol kinase family.,similarity:Contains 1 DAGKc domain.,similarity:Contains 2 EF-hand domains.,similarity:Contains 2 phorbol-ester/DAG-type zinc fingers.,


研究领域 (Research Area)

Glycerolipid metabolism;Glycerophospholipid metabolism;Phosphatidylinositol signaling system;

图片 (Image Data)



Western blot analysis of the lysates from HUVEC cells using DGKB antibody.



Western Blot analysis of HuvEC cells using DGK- β Polyclonal Antibody

注意事项 (Note)

For research use only .

产品名称: DGK- β Rabbit Polyclonal Antibody
产品货号: APRab09947

