

产品名称: DHRS2 Rabbit Polyclonal Antibody
产品货号: APRab09970



产品概述 (Summary)

产品名称 (Production Name)	DHRS2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DHRS2
别名 (Alternative Names)	DHRS2; Dehydrogenase/reductase SDR family member 2; Dicarbonyl reductase HEP27; Protein D
基因 ID (Gene ID)	10202.0
蛋白 ID (SwissProt ID)	Q13268.The antiserum was produced against synthesized peptide derived from human DHRS2. AA range:111-160

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	27kDa

产品名称: DHRS2 Rabbit Polyclonal Antibody
产品货号: APRab09970




研究背景 (Background)


function: May inhibit cell replication either by catalyzing the oxidation of estrogen and androgen or by converting cortisone in cortisol., similarity: Belongs to the short-chain dehydrogenases/reductases (SDR) family., function: May inhibit cell replication either by catalyzing the oxidation of estrogen and androgen or by converting cortisone in cortisol., similarity: Belongs to the short-chain dehydrogenases/reductases (SDR) family.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from COLO cells, using DHRS2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HeLa cells using DHRS2 antibody.

产品名称: DHRS2 Rabbit Polyclonal Antibody
产品货号: AP Rab09970



Western Blot analysis of various cells using DHRS2 Polyclonal Antibody

注意事项 (Note)

For research use only .