

产品名称: DLEU7 Rabbit Polyclonal Antibody
产品货号: APRab10014



产品概述 (Summary)

产品名称 (Production Name)	DLEU7 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DLEU7 LEU7
别名 (Alternative Names)	Leukemia-associated protein 7 (Deleted in lymphocytic leukemia 7)
基因 ID (Gene ID)	220107.0
蛋白 ID (SwissProt ID)	Q6UYE1. The antiserum was produced against synthesized peptide derived from the Internal region of human DLEU7. AA range: 51-100

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:50-1:200, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: DLEU7 Rabbit Polyclonal Antibody
产品货号: APRab10014



alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and kidney.,tissue specificity:Ubiquitous.,alternative products:Additional isoforms seem to exist,developmental stage:Expressed in the fetal brain, liver and kidney.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

注意事项 (Note)

For research use only .