

产品名称: DQX1 Rabbit Polyclonal Antibody
产品货号: APRab10148



产品概述 (Summary)

产品名称 (Production Name)	DQX1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DQX1
别名 (Alternative Names)	DQX1; ATP-dependent RNA helicase DQX1; DEAQ box polypeptide 1
基因 ID (Gene ID)	165545.0
蛋白 ID (SwissProt ID)	Q8TE96.The antiserum was produced against synthesized peptide derived from human DQX1. AA range:551-600

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	79kDa

研究背景 (Background)


产品名称: DQX1 Rabbit Polyclonal Antibody
产品货号: APRab10148




similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,similarity:Contains 1 helicase ATP-binding domain.,similarity:Contains 1 helicase C-terminal domain.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from COS7 cells, using DQX1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western Blot analysis of various cells using DQX1 Polyclonal Antibody cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Invent biotech, MN, USA) .

注意事项 (Note)

For research use only .