

产品名称: Dynein LC 2B Rabbit Polyclonal Antibody
产品货号: APRab10231



产品概述 (Summary)

产品名称 (Production Name)	Dynein LC 2B Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DYNLRB2
别名 (Alternative Names)	DYNLRB2; DNCL2B; DNLC2B; ROBLD2; Dynein light chain roadblock-type 2; Dynein light chain 2B; cytoplasmic; Roadblock domain-containing protein 2
基因 ID (Gene ID)	83657.0
蛋白 ID (SwissProt ID)	Q8TF09. The antiserum was produced against synthesized peptide derived from human DYNLRB2. AA range: 21-70

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	

产品名称: Dynein LC 2B Rabbit Polyclonal Antibody
产品货号: APRab10231




研究背景 (Background)

function: May be involved in assembly and motor function of dynein, which plays a central role in cell division and intracellular transport., miscellaneous: Expression is significantly down-regulated in hepatocellular carcinoma (HCC) patients., similarity: Belongs to the GAMAD family., tissue specificity: High expression in heart, brain, placenta, skeletal muscle, prostate and small intestine; moderate in kidney, pancreas, spleen, testis, ovary and colon; low in lung, liver, thymus and leukocyte., function: May be involved in assembly and motor function of dynein, which plays a central role in cell division and intracellular transport., miscellaneous: Expression is significantly down-regulated in hepatocellular carcinoma (HCC) patients., similarity: Belongs to the GAMAD family., tissue specificity: High expression in heart, brain, placenta, skeletal muscle, prostate and small intestine; moderate in kidney, pancreas, spleen, testis, ovary and colon; low in lung, liver, thymus and leukocyte.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human testis tissue, using DYNLRB2 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .