

产品名称: DZIP3 Rabbit Polyclonal Antibody
产品货号: APRab10248



产品概述 (Summary)

产品名称 (Production Name)	DZIP3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	DZIP3
别名 (Alternative Names)	DZIP3; KIAA0675; E3 ubiquitin-protein ligase DZIP3; DAZ-interacting protein 3; RNA-binding ubiquitin ligase of 138 kDa; hRUL138
基因 ID (Gene ID)	9666.0
蛋白 ID (SwissProt ID)	Q86Y13.The antiserum was produced against synthesized peptide derived from human DZIP3. AA range:681-730

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	140kDa

产品名称: DZIP3 Rabbit Polyclonal Antibody
产品货号: AP Rab10248




研究背景 (Background)

function: E3 Ubiquitin ligase proteins mediate ubiquitination and subsequent proteasomal degradation of target proteins. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Able to specifically bind RNA., pathway: Protein modification; protein ubiquitination., sequence caution: Contaminating sequence. Potential poly-A sequence., similarity: Contains 1 RING-type zinc finger., subunit: Interacts with DAZ proteins., tissue specificity: Widely expressed at low level. Highly expressed in skeletal muscle, kidney and heart. Expressed at low level in placenta, lung, brain, liver and pancreas., function: E3 Ubiquitin ligase proteins mediate ubiquitination and subsequent proteasomal degradation of target proteins. E3 ubiquitin ligases accept ubiquitin from an E2 ubiquitin-conjugating enzyme in the form of a thioester and then directly transfers the ubiquitin to targeted substrates. Able to specifically bind RNA., pathway: Protein modification; protein ubiquitination., sequence caution: Contaminating sequence. Potential poly-A sequence., similarity: Contains 1 RING-type zinc finger., subunit: Interacts with DAZ proteins., tissue specificity: Widely expressed at low level. Highly expressed in skeletal muscle, kidney and heart. Expressed at low level in placenta, lung, brain, liver and pancreas.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of DZIP3 Antibody. The lane on the right is blocked with the DZIP3 peptide.

注意事项 (Note)

For research use only .