

产品名称: Ebi3 Rabbit Polyclonal Antibody
产品货号: APRab10273



产品概述 (Summary)

产品名称 (Production Name)	Ebi3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EBI3
别名 (Alternative Names)	EBI3; IL27B; Interleukin-27 subunit beta; IL-27 subunit beta; IL-27B; Epstein-Barr virus-induced gene 3 protein; EBV-induced gene 3 protein
基因 ID (Gene ID)	10148.0
蛋白 ID (SwissProt ID)	Q14213.The antiserum was produced against synthesized peptide derived from the Internal region of human EBI3. AA range:101-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	25kDa


产品名称: Ebi3 Rabbit Polyclonal Antibody
产品货号: AP Rab10273

研究背景 (Background)

This gene was identified by its induced expression in B lymphocytes in response Epstein-Barr virus infection. It encodes a secreted glycoprotein belonging to the hematopoietin receptor family, and heterodimerizes with a 28 kDa protein to form interleukin 27 (IL-27). IL-27 regulates T cell and inflammatory responses, in part by activating the Jak/STAT pathway of CD4⁺ T cells. [provided by RefSeq, Sep 2008],function:Cytokine with pro- and anti-inflammatory properties, that can regulate T helper cell development, suppress T-cell proliferation, stimulate cytotoxic T cell activity, induce isotype switching in B-cells, and that has diverse effects on innate immune cells. Among its target cells are CD4 T helper cells which can differentiate in type 1 effector cells (TH1), type 2 effector cells (TH2) and IL17 producing helper T-cells (TH17). It drives rapid clonal expansion of naive but not memory CD4 T-cells. It also strongly synergizes with IL-12 to trigger interferon-gamma/IFN-gamma production of naive CD4 T-cells, binds to the cytokine receptor WSX-1/TCCR. Another important role of IL27 is its antitumor activity as well as its antiangiogenic activity with activation of production of antiangiogenic chemokines.,induction:By EBV.,online information:Interleukin-27 entry,similarity:Belongs to the type I cytokine receptor family. Type 3 subfamily.,similarity:Contains 2 fibronectin type-III domains.,subunit:Heterodimer with IL27. IL-27 is an heterodimer composed IL27 and EBI3. EBI3 is also a component of the IL-12 heterodimer. Interacts with SQSTM1.,

研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of HBE cells using Ebi3 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

注意事项 (Note)

For research use only .