

产品名称: EBP2 Rabbit Polyclonal Antibody
产品货号: APRab10276



产品概述 (Summary)

产品名称 (Production Name)	EBP2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EBNA1BP2 EBP2
别名 (Alternative Names)	
基因 ID (Gene ID)	10969.0
蛋白 ID (SwissProt ID)	Q99848.Synthesized peptide derived from human protein . at AA range: 140-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	33kDa

研究背景 (Background)

function:Required for the processing of the 27S pre-rRNA.,similarity:Belongs to the EBP2 family.,subunit:Specifically


产品名称: EBP2 Rabbit Polyclonal Antibody
产品货号: APRab10276



interacts with EBV EBNA1. The EBNA1-EBP2 interaction is important for the stable segregation of EBV episomes during cell division.,tissue specificity:Ubiquitous.,function:Required for the processing of the 27S pre-rRNA.,similarity:Belongs to the EBP2 family.,subunit:Specifically interacts with EBV EBNA1. The EBNA1-EBP2 interaction is important for the stable segregation of EBV episomes during cell division.,tissue specificity:Ubiquitous.,

研究领域 (Research Area)

图片 (Image Data)



Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4°over night

注意事项 (Note)

For research use only .