

产品名称: ECP Rabbit Polyclonal Antibody
产品货号: APRab10285



产品概述 (Summary)

产品名称 (Production Name)	ECP Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	RNASE3
别名 (Alternative Names)	RNASE3; ECP; RNS3; Eosinophil cationic protein; ECP; Ribonuclease 3; RNase 3
基因 ID (Gene ID)	6037.0
蛋白 ID (SwissProt ID)	P12724.The antiserum was produced against synthesized peptide derived from human ECP. AA range:102-151

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	18kDa

研究背景 (Background)

产品名称: ECP Rabbit Polyclonal Antibody
产品货号: APRab10285




Human, Mouse, Rat, Rabbit, Fish, Chicken, Guinea pig, Sheep, Bovine


研究领域 (Research Area)

Asthma;

图片 (Image Data)



Western blot analysis of lysate from 293 cells, using ECP antibody.



Western Blot analysis of various cells using ECP Polyclonal Antibody

注意事项 (Note)

For research use only .