

产品名称: EF-1 γ Rabbit Polyclonal Antibody
产品货号: APRab10313



产品概述 (Summary)

产品名称 (Production Name)	EF-1 γ Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1G
别名 (Alternative Names)	EEF1G; EF1G; Elongation factor 1-gamma; EF-1-gamma; eEF-1B gamma
基因 ID (Gene ID)	1937.0
蛋白 ID (SwissProt ID)	P26641.The antiserum was produced against synthesized peptide derived from human EEF1G. AA range:101-150

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	55kDa

研究背景 (Background)



产品名称: EF-1 γ Rabbit Polyclonal Antibody
产品货号: APRab10313




eukaryotic translation elongation factor 1 gamma (EEF1G) Homo sapiens This gene encodes a subunit of the elongation factor-1 complex, which is responsible for the enzymatic delivery of aminoacyl tRNAs to the ribosome. This subunit contains an N-terminal glutathione transferase domain, which may be involved in regulating the assembly of multisubunit complexes containing this elongation factor and aminoacyl-tRNA synthetases. [provided by RefSeq, Jul 2008], function: Probably plays a role in anchoring the complex to other cellular components., similarity: Contains 1 EF-1-gamma C-terminal domain., similarity: Contains 1 GST C-terminal domain., similarity: Contains 1 GST N-terminal domain., subunit: EF-1 is composed of four subunits: alpha, beta, delta, and gamma., tissue specificity: Highly expressed in pancreatic tumor tissue and to a lesser extent in normal kidney, intestine, pancreas, stomach, lung, brain, spleen and liver.,

研究领域 (Research Area)

图片 (Image Data)



产品名称: EF-1 γ Rabbit Polyclonal Antibody
产品货号: APRab10313



Western Blot analysis of various cells using EF-1 γ Polyclonal Antibody diluted at 1: 1000

注意事项 (Note)

For research use only .