

产品名称: EF-1 β Rabbit Polyclonal Antibody
产品货号: APRab10317



产品概述 (Summary)

产品名称 (Production Name)	EF-1 β Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Rat

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EEF1B2
别名 (Alternative Names)	EEF1B2; EEF1B; EF1B; Elongation factor 1-beta; EF-1-beta
基因 ID (Gene ID)	1933.0
蛋白 ID (SwissProt ID)	P24534.The antiserum was produced against synthesized peptide derived from human EF-1beta. AA range:83-132

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	24kDa

研究背景 (Background)


产品名称: EF-1 β Rabbit Polyclonal Antibody
产品货号: APRab10317



This gene encodes a translation elongation factor. The protein is a guanine nucleotide exchange factor involved in the transfer of aminoacylated tRNAs to the ribosome. Alternative splicing results in three transcript variants which differ only in the 5' UTR. [provided by RefSeq, Jul 2008],function:EF-1-beta and EF-1-delta stimulate the exchange of GDP bound to EF-1-alpha to GTP.,PTM:Phosphorylation affects the GDP/GTP exchange rate.,similarity:Belongs to the EF-1-beta/EF-1-delta family.,similarity:Contains 1 GST C-terminal domain.,subunit:EF-1 is composed of 4 subunits: alpha, beta, delta, and gamma.,

研究领域 (Research Area)

图片 (Image Data)



注意事项 (Note)

For research use only .