

产品名称: ELL2 Rabbit Polyclonal Antibody
产品货号: APRab10407



产品概述 (Summary)

产品名称 (Production Name)	ELL2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	IHC, ICC, IF, ELISA
种属反应性 (Reactivity)	Human, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ELL2
别名 (Alternative Names)	ELL2; RNA polymerase II elongation factor ELL2
基因 ID (Gene ID)	22936.0
蛋白 ID (SwissProt ID)	O00472. The antiserum was produced against synthesized peptide derived from human ELL2. AA range: 448-497

产品应用 (Application)

稀释比 (Dilution Ratio)	IHC 1:100-1:300, ICC/IF 1:50-1:200, ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	

研究背景 (Background)


产品名称: ELL2 Rabbit Polyclonal Antibody
产品货号: APRab10407



function:Elongation factor that can increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.,similarity:Belongs to the ELL/occludin family.,subunit:Interacts with EAF1 and EAF2.,function:Elongation factor that can increase the catalytic rate of RNA polymerase II transcription by suppressing transient pausing by the polymerase at multiple sites along the DNA.,similarity:Belongs to the ELL/occludin family.,subunit:Interacts with EAF1 and EAF2.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma tissue, using ELL2 Antibody. The picture on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .