

产品名称: EPG5 Rabbit Polyclonal Antibody
产品货号: APRab10509



产品概述 (Summary)

产品名称 (Production Name)	EPG5 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF
种属反应性 (Reactivity)	Human,Rat,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EPG5 KIAA1632
别名 (Alternative Names)	Ectopic P granules protein 5 homolog
基因 ID (Gene ID)	57724.0
蛋白 ID (SwissProt ID)	Q9HCE0.Synthetic peptide from human protein at AA range: 141-190

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:1000-1:2000,IHC 1:200-1:500,ICC/IF 1:50-1:200
蛋白分子量 (Molecular Weight)	

研究背景 (Background)

This gene encodes a large coiled coil domain-containing protein that functions in autophagy during starvation conditions.


产品名称: EPG5 Rabbit Polyclonal Antibody
产品货号: AP Rab10509




Mutations in this gene cause Vici syndrome. [provided by RefSeq, Aug 2015],sequence caution:Translated as Trp.,similarity:Belongs to the UPF0493 family.,

研究领域 (Research Area)

图片 (Image Data)



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-liver, antibody was diluted at 1:200

注意事项 (Note)

For research use only .