

产品名称: Eps8L2 Rabbit Polyclonal Antibody
产品货号: APRab10549



产品概述 (Summary)

产品名称 (Production Name)	Eps8L2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	EPS8L2 EPS8L2; EPS8R2; PP13181; Epidermal growth factor receptor kinase substrate
别名 (Alternative Names)	8-like protein 2; EPS8-like protein 2; Epidermal growth factor receptor pathway substrate 8-related protein 2; EPS8-related protein 2
基因 ID (Gene ID)	64787.0
蛋白 ID (SwissProt ID)	Q9H6S3.The antiserum was produced against synthesized peptide derived from human ES8L2. AA range:1-50

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:20000-1:40000
蛋白分子量 (Molecular Weight)	81kDa

产品名称: Eps8L2 Rabbit Polyclonal Antibody
产品货号: APRab10549





研究背景 (Background)

This gene encodes a member of the EPS8 gene family. The encoded protein, like other members of the family, is thought to link growth factor stimulation to actin organization, generating functional redundancy in the pathways that regulate actin cytoskeletal remodeling. [provided by RefSeq, Dec 2008],function:Stimulates guanine exchange activity of SOS1. May play a role in membrane ruffling and remodeling of the actin cytoskeleton.,similarity:Belongs to the EPS8 family.,similarity:Contains 1 PID domain.,similarity:Contains 1 SH3 domain.,subunit:Interacts with ABI1. Part of a complex that contains SOS1, ABI1 and EPS8L2. Associates with F-actin.,tissue specificity:Detected in fibroblasts and placenta.,


研究领域 (Research Area)

图片 (Image Data)



Western Blot analysis of various cells using Eps8L2 Polyclonal Antibody

产品名称: Eps8L2 Rabbit Polyclonal Antibody
产品货号: APRab10549



Western Blot analysis of HuvEc cells using Eps8L2 Polyclonal Antibody

注意事项 (Note)

For research use only .