

产品名称: ERGIC-3 Rabbit Polyclonal Antibody
产品货号: APRab10590



产品概述 (Summary)

产品名称 (Production Name)	ERGIC-3 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,IHC,ICC/IF,ELISA
种属反应性 (Reactivity)	Human,Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	ERGIC3 ERGIC3; C20orf47; ERV46; SDBCAG84; CGI-54; Endoplasmic reticulum-Golgi intermediate compartment protein 3; Serologically defined breast cancer antigen NY-BR-84
别名 (Alternative Names)	
基因 ID (Gene ID)	51614.0
蛋白 ID (SwissProt ID)	Q9Y282.The antiserum was produced against synthesized peptide derived from human ERG13. AA range:321-370

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:200-1:1000,ELISA 1:10000-1:20000
蛋白分子量 (Molecular Weight)	43kDa

产品名称: ERGIC-3 Rabbit Polyclonal Antibody
产品货号: APRab10590




研究背景 (Background)


function: Possible role in transport between endoplasmic reticulum and Golgi., similarity: Belongs to the ERGIC family., subcellular location: Cycles between the endoplasmic reticulum and the Golgi., subunit: Interacts with ERGIC1/ERGIC32., function: Possible role in transport between endoplasmic reticulum and Golgi., similarity: Belongs to the ERGIC family., subcellular location: Cycles between the endoplasmic reticulum and the Golgi., subunit: Interacts with ERGIC1/ERGIC32.,

研究领域 (Research Area)

图片 (Image Data)



Immunofluorescence analysis of HepG2 cells, using ERGIC3 Antibody. The picture on the right is blocked with the synthesized



Western blot analysis of lysates from LOVO and NIH/3T3 cells, using ERGIC3 Antibody. The lane on the right is blocked with the synthesized peptide.

注意事项 (Note)

For research use only .

产品名称: ERGIC-3 Rabbit Polyclonal Antibody
产品货号: APRab10590

