

**产品名称: FAF1 Rabbit Polyclonal Antibody**  
**产品货号: APRab10799**



## 产品概述 (Summary)

|                        |                                 |
|------------------------|---------------------------------|
| 产品名称 (Production Name) | FAF1 Rabbit Polyclonal Antibody |
| 描述 (Description)       | Rabbit polyclonal Antibody      |
| 宿主 (Host)              | Rabbit                          |
| 应用 (Application)       | WB,IHC,ICC/IF,ELISA             |
| 种属反应性 (Reactivity)     | Human,Mouse,Rat                 |

## 产品性能 (Performance)

|                     |                                                                                                   |
|---------------------|---------------------------------------------------------------------------------------------------|
| 偶联物 (Conjugation)   | Unconjugated                                                                                      |
| 修饰 (Modification)   | Unmodified                                                                                        |
| 同种型 (Isotype)       | IgG                                                                                               |
| 克隆 (Clonality)      | Polyclonal                                                                                        |
| 形式 (Form)           | Liquid                                                                                            |
| 存放说明 (Storage)      | Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.          |
| 储存溶液 (Buffer)       | Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N. |
| 纯化方式 (Purification) | Affinity purification                                                                             |

## 免疫原信息 (Immunogen)

|                        |                                                                                                                                  |
|------------------------|----------------------------------------------------------------------------------------------------------------------------------|
| 基因名 (Gene Name)        | FAF1                                                                                                                             |
| 别名 (Alternative Names) | FAF1; UBXD12; UBXN3A; CGI-03; FAS-associated factor 1; hFAF1; UBX domain-containing protein 12; UBX domain-containing protein 3A |
| 基因 ID (Gene ID)        | 11124.0                                                                                                                          |
| 蛋白 ID (SwissProt ID)   | Q9UNN5.The antiserum was produced against synthesized peptide derived from the Internal region of human FAF1. AA range:461-510   |

## 产品应用 (Application)

|                          |                                                                         |
|--------------------------|-------------------------------------------------------------------------|
| 稀释比 (Dilution Ratio)     | WB 1:500-1:2000,IHC 1:100-1:300,ICC/IF 1:50-1:200,ELISA 1:10000-1:20000 |
| 蛋白分子量 (Molecular Weight) | 75kDa                                                                   |

产品名称: FAF1 Rabbit Polyclonal Antibody  
产品货号: APRab10799




## 研究背景 (Background)


Interaction of Fas ligand (TNFSF6) with the FAS antigen (TNFRSF6) mediates programmed cell death, also called apoptosis, in a number of organ systems. The protein encoded by this gene binds to FAS antigen and can initiate apoptosis or enhance apoptosis initiated through FAS antigen. Initiation of apoptosis by the protein encoded by this gene requires a ubiquitin-like domain but not the FAS-binding domain. [provided by RefSeq, Jul 2008],function:Potentiates but cannot initiate FAS-induced apoptosis.,similarity:Contains 1 UBX domain.,subunit:Specifically interacts with the cytoplasmic domain of FAS.,tissue specificity:Most abundant in testis, slightly less abundant in skeletal muscle and heart, followed by prostate, thymus, ovary, small intestine, and colon. Not detected in the peripheral blood leukocytes.,

## 研究领域 (Research Area)

## 图片 (Image Data)




Western blot analysis of lysate from K562 cells, using FAF1 Antibody.



Western Blot analysis of K562 cells using FAF1 Polyclonal Antibody.. Secondary antibody was diluted at 1:20000

**产品名称: FAF1 Rabbit Polyclonal Antibody**  
**产品货号: APRab10799**



Immunohistochemical analysis of paraffin-embedded mouse-brain, antibody was diluted at 1:100

### **注意事项 (Note)**

For research use only .