

产品名称: FASTKD1 Rabbit Polyclonal Antibody
产品货号: APRab10842



产品概述 (Summary)

产品名称 (Production Name)	FASTKD1 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB,ELISA
种属反应性 (Reactivity)	Human,Mouse,Monkey

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FASTKD1
别名 (Alternative Names)	FASTKD1; KIAA1800; FAST kinase domain-containing protein 1
基因 ID (Gene ID)	79675.0
蛋白 ID (SwissProt ID)	Q53R41.The antiserum was produced against synthesized peptide derived from human FAKD1. AA range:561-610

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000,ELISA 1:5000-1:20000
蛋白分子量 (Molecular Weight)	97kDa

研究背景 (Background)


产品名称: FASTKD1 Rabbit Polyclonal Antibody
产品货号: APRab10842




caution:It is uncertain whether Met-1 or Met-15 is the initiator.,similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,caution:It is uncertain whether Met-1 or Met-15 is the initiator.,similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,

研究领域 (Research Area)


图片 (Image Data)



Western blot analysis of lysates from COS7 cells, using FASTKD1 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from HT-29 cells using FASTKD1 antibody.



产品名称: FASTKD1 Rabbit Polyclonal Antibody
产品货号: APRab10842



Western Blot analysis of various cells using FASTKD1 Polyclonal Antibody

注意事项 (Note)

For research use only .