

产品名称: FASTKD2 Rabbit Polyclonal Antibody
产品货号: APRab10843



产品概述 (Summary)

产品名称 (Production Name)	FASTKD2 Rabbit Polyclonal Antibody
描述 (Description)	Rabbit polyclonal Antibody
宿主 (Host)	Rabbit
应用 (Application)	WB, ICC/IF, ELISA
种属反应性 (Reactivity)	Human, Rat, Mouse

产品性能 (Performance)

偶联物 (Conjugation)	Unconjugated
修饰 (Modification)	Unmodified
同种型 (Isotype)	IgG
克隆 (Clonality)	Polyclonal
形式 (Form)	Liquid
存放说明 (Storage)	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
储存溶液 (Buffer)	Liquid in PBS containing 50% glycerol, 0.5% protective protein and 0.02% New type preservative N.
纯化方式 (Purification)	Affinity purification

免疫原信息 (Immunogen)

基因名 (Gene Name)	FASTKD2
别名 (Alternative Names)	FASTKD2; KIAA0971; FAST kinase domain-containing protein 2
基因 ID (Gene ID)	22868.0
蛋白 ID (SwissProt ID)	Q9NYY8. The antiserum was produced against synthesized peptide derived from human FAKD2. AA range: 171-220

产品应用 (Application)

稀释比 (Dilution Ratio)	WB 1:500-1:2000, ICC/IF 1:50-1:200, ELISA 1:5000-1:10000
蛋白分子量 (Molecular Weight)	81kDa

研究背景 (Background)


产品名称: FASTKD2 Rabbit Polyclonal Antibody
产品货号: AP Rab10843




This gene encodes a protein that is localized in the mitochondrial inner compartment and that may play a role in mitochondrial apoptosis. Nonsense mutations have been reported to result in cytochrome c oxidase deficiency. [provided by RefSeq, Oct 2008],caution:It is uncertain whether Met-1 or Met-17 is the initiator.,similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,

研究领域 (Research Area)

图片 (Image Data)




Western blot analysis of lysates from HUVEC, HepG2, and COLO cells, using FAKD2 Antibody. The lane on the right is blocked with the synthesized peptide.



Western blot analysis of the lysates from COLO205 cells using FAKD2 antibody.

产品名称: FASTKD2 Rabbit Polyclonal Antibody
产品货号: APRab10843



Western Blot analysis of various cells using FASTKD2 Polyclonal Antibody

注意事项 (Note)

For research use only .